


# Mikki





Padding that grows with a child



Trunk and pelvis guide as well as adjustable seat depth



Adjustable back height

# Mikki



# The Best Posture

# Options

## **Trunk supports**

The height and width of the thoracic supports can be adjusted and offer additional postural support to the trunk.



## **Adduction guide**

The width and angle of the gluteal supports and adduction guides can be set. It is also possible to combine different sizes with each other.



## **Abductor**

A swing-away or detachable abductor is also available as an option. This feature reduces the adduction of the legs and holds the hips and thigh in a neutral position.



# Mikki



## Mikki at a glance

- Active seating shell
- Mikki is attached to the wheelchair with belts
- Available in three sizes; adjustable parameters, see sizing chart
- Pre-shaped padding elements facilitate individual adaptation
- Shifting lumbar pad
- All pads are individually adjustable
- Abduction from 0° to 4° on each side can be set
- Option: swing-away and detachable abductor
- Different-sized seat and back pieces can be combined
- Cover is black, incontinence-proof, removable and washable at 30° C
- Safety systems



## Technical data in cm and kg



	Size 1	Size 2	Size 3
Seat width	18-24	24-31	30-40
Seat depth	18-26	24-32	31-41
Height of backrest	19-26	25-32	31-42
Axle height	16-22	21-29	28-40
Thoracic width	17-23	23-30	28-39